



62179

## Aprotinin ex. Bovine Lung, 5400U/mg

Part B

CAS : 9087-70-1

### Specifications

Appearance (Colour)	White
Appearance (Form)	Powder
Activity	min. 5400 units/mg powder
Unit Definition	One trypsin inhibitor unit (TIU) will decrease the activity of two trypsin units by 50% where one trypsin unit will hydrolyze 1.0 micromole of N <sup>?</sup> -Benzoyl-DL-Arginine-p-Nitroaniline (BAPNA) per minute at pH 7.8 at 25 °C.
Enzyme solution for activity determination	Prepared in 1 mM hydrochloric acid solution.

### General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Mg	35079099 (GST 18%)
100 Mg	35079099 (GST 18%)
25 Mg	35079099 (GST 18%)
5 Mg	35079099 (GST 18%)

### Available Packages

5 Mg

---

10 Mg

---

25 Mg

---

100 Mg

---

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com