



62205

Uridine extrapure, 99%

Part B

CAS : 58-96-8

Molecular Formula : C₉H₁₂N₂O₆

Molecular Weight : 244.20

Specifications

Appearance (Colour)	White to off - white
Appearance (Form)	Crystalline powder
Assay (UV)	min.99%
Melting Point	166-167°C
Loss on drying	max.0.5%
Heavy Metals (Pb)	max. 0.001%
Spectral Ratio (A250nm/A260nm)	~0.74
Spectral Ratio (A280nm/A260nm)	~0.35

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29349990 (GST 18%)
25 Gms	29349990 (GST 18%)
5 Gms	29349990 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
25 Gms	Plastic Bottle

Available Packages

5 Gms
25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com