



62399

## Sodium Selenate Anhydrous extrapure AR, 99%

Part A

CAS : 13410-01-0

Molecular Formula : Na<sub>2</sub>SeO<sub>4</sub>

Molecular Weight : 188.94

### Specifications

Appearance (Colour)	White to Off-white
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear to slight opalescent
Assay (Iodometric)	min.99%
Loss on drying	max.1%
Chloride (Cl)	max.0.005%
Sulphate (SO <sub>4</sub> )	max.0.01%
Iron (Fe)	max.0.001%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2630 , IMCO Class No.:6.1 , Packing Group:I
HSN Code	
100 Gms	28429090 (GST 18%)
25 Gms	28429090 (GST 18%)
500 Gms	28429090 (GST 18%)

### Available Packages

25 Gms  
100 Gms  
500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com