



62603

## Lead Metal Powder extrapure AR, 99.5%

Part B

CAS : 7439-92-1

Molecular Formula : Pb

Molecular Weight : 207.20

### Specifications

Appearance (Colour)	Greyish black
Appearance (Form)	Fine powder
Solubility (Turbidity) 5% solution in 5N HNO <sub>3</sub>	Clear
Assay (EDTA Titration)	min. 99.5%
Iron (Fe)	max. 0.001%
Copper (Cu)	max. 0.002%
Arsenic (As)	max. 0.0001%
Silver (Ag)	max. 0.002%
Tin (Sn)	max. 0.001%
Antimony (Sb)	max. 0.005%
Nickel (Ni)	max. 0.0005%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3077 , IMCO Class No.:9 , Packing Group:III
HSN Code	
500 Gms	78042000 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

### Available Packages

500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com