



63338

## Silicon Nitride Nanopowder, 20nm, 99%

Part C

CAS : 12033-89-5

Molecular Formula : Si<sub>3</sub>N<sub>4</sub>

Molecular Weight : 140.28

### Specifications

Appearance (Colour)	Off white to beige to grey
Appearance (Form)	Powder
Assay	min. 99%
APS	20nm

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
1 Gms	28500020 (GST 18%)
5 Gms	28500020 (GST 18%)
Type of Packing	
5 Gms	Plastic Bottle

### Available Packages

5 Gms  
25 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com