



63390

Yttrium (III) Chloride Hexahydrate extrapure AR, 99.9%

Part B

CAS : 10025-94-2

Molecular Formula : $YCl_3 \cdot 6H_2O$

Molecular Weight : 303.36

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline compound
Assay (EDTA Titration)	min. 99.9%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	28469010 (GST 18%)
25 Gms	28469010 (GST 18%)
5 Gms	28469010 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

Available Packages

5 Gms
25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com