



63504

Salicylaldehyde extrapure AR, ExiPlus, Multi-Compendial, 99.5%

Part A

CAS : 90-02-8

Molecular Formula : C7H6O2

Molecular Weight : 122.12

Specifications

ExiPlus	meets compendial specs of Ph.Eur
Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay (GC)	min. 99.5%
Density (g/ml) @ 20°C	1.164-1.167
Refractive Index (20°C)	1.573-1.574
Boiling Range	194-197°C
Melting Point	2-5°C
Acidity	max. 0.1% (C7H6O3)
Water (KF)	max. 0.05%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	29122990 (GST 18%)
250 ml	29122990 (GST 18%)
500 ml	29122990 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

Available Packages

100 ml

250 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com