



63661

Indole-3-Acetyl-L-Phenylalanine extrapure, 98%

Part B

CAS : 57105-50-7

Molecular Formula : C₁₉H₁₈N₂O₃

Molecular Weight : 322.37

General Information

Storage 2 to 8°C (Refrigerate)

HSN Code

1 Gms

98020000 (GST 18%)

5 Gms

98020000 (GST 18%)

Available Packages

1 Gms

5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com