



63703

Sodium Taurocholate Hydrate (STC Hydrate) extrapure, 98%

Part B

CAS : 345909-26-4

Molecular Formula : C₂₆H₄₅NNaO₇S.xH₂O

Molecular Weight : 537.68 (anhy)

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 5% aq. solution	Clear
Solubility (Colour) 5% aq. solution	Colourless
Assay (On dry basis)	min. 98%
Water (KF)	max. 5.0%

Other Information

Application

Bile Salt suitable for preparation of media and simulated gastric and intestinal fluids for Pharmaceutical Quality Control applications.

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29242990 (GST 18%)
100 Gms	29242990 (GST 18%)
25 Gms	29242990 (GST 18%)
5 Gms	29242990 (GST 18%)
Type of Packing	
1 Gms	Glass Bottle
25 Gms	Glass Bottle

Available Packages

1 Gms

5 Gms

25 Gms

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
