



63797

Xylene Cyanol FF

Part A

CAS : 2650-17-1

Molecular Formula : C₂₅H₂₇N₂NaO₆S₂

Molecular Weight : 538.6

Specifications

Appearance (Colour)	Green
Appearance (Form)	Crystalline powder
Loss on drying	max. 8%
Absorbance (A) of 1% aq. solution in a 1cm cell v/s H ₂ O	
@614nm	1000 - 1400
Molar Extinction coefficient (?) in Water at 613 - 616 nm	min. 65000(On dried basis)

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	32041248 (GST 18%)
25 Gms	32041248 (GST 18%)
5 Gms	32041248 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle

Available Packages

5 Gms

25 Gms

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com