



63952

Butyl Propionate extrapure, 99%

Part A
CAS : 590-01-2
Molecular Formula : C7H14O2
Molecular Weight : 130.19

Specifications	
Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 25
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	0.875 - 0.877
Refractive Index (20°C)	1.399 - 1.401
Boiling Range	142 - 145°C
Water (KF)	max. 0.1%

General Information	
Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No: 1914 , IMCO Class No.: 3 , Packing Group:III
HSN Code	
500 ml	29155000 (GST 18%)
2500 ml	29155000 (GST 18%)

Available Packages	
500 ml	
2500 ml	

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.
608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com