



64187

Magnesium Sulphate Dried extrapure AR, ExiPlus, Multi-Compendial, 62-70% MgSO₄, 99%

Part A

CAS : 7487-88-9

Molecular Formula : MgSO₄.xH₂O

Molecular Weight : 120.37 (anhy)

Specifications

ExiPlus	meets compendial specs of USP, BP, Ph.Eur
Appearance (Colour)	White
Appearance (Form)	Crystalline compound
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (as is)	62 - 70%
Assay (On dry basis)	min. 99%
Loss on drying	30 -35%
Chloride (Cl)	max. 0.001%
Iron (Fe)	max. 0.001%
Heavy Metals (Pb)	max. 0.001%
Nitrate (NO ₃)	max. 0.004%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
2.5 Kg	28332940 (GST 18%)
500 Gms	28332940 (GST 18%)
Type of Packing	
2.5 Kg	Plastic Jar
500 Gms	Plastic Bottle

Available Packages

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com