



64674

Tetraethylammonium Hydrogen Sulphate
(TEAHS) extrapure, 99%

Part A

CAS : 16873-13-5

Molecular Formula : (C₂H₅)₄N(Cl)

Molecular Weight : 165.70

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (T)	min. 99%
Absorbance (A) of 0.005 M aq. solution in a 1 cm cell v/s H ₂ O	
@200nm	max. 0.15
@220nm	max. 0.05
@250 nm	max. 0.01

Other Information

Application	Suitable as an Ion Pairing Reagent & for HPLC applications.
-------------	---

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	29239000 (GST 18%)
25 Gms	29239000 (GST 18%)

Available Packages

10 Gms
25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com