



64839

Bilirubin extrapure AR, 99%

Part B

CAS : 635-65-4

Molecular Formula : C₃₃H₃₆N₄O₆

Molecular Weight : 584.68

Specifications

Appearance (Colour)	Yellowish orange
Appearance (Form)	Crystalline powder
Assay	min. 99%
Absorbance (A) of 1% solution in Chloroform in a 1cm cell @454nm	min.1030
Loss on drying	max. 0.1%
Sulphated Ash	max. 0.1%

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29337990 (GST 18%)
100 Mg	29337990 (GST 18%)
Type of Packing	
100 Mg	Glass Bottle

Available Packages

100 Mg

1 Gms

5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com