



64857

Iron Metal Powder (Electrolytic), 99.5%, ~200 mesh

Part A

CAS : 7439-89-6

Molecular Formula : Fe

Molecular Weight : 55.85

Specifications

Appearance (Colour)	Blackish to grey
Appearance (Form)	Fine powder
Solubility (Turbidity) 5% solution in 5N HCl	Clear
Assay (Redox Titration)	min. 99.5%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3089 , IMCO Class No.:4.1 , Packing Group:III
HSN Code	
1 Kg	72052910 (GST 18%)
5 Kg	72052910 (GST 18%)
500 Gms	72052910 (GST 18%)
Type of Packing	
1 Kg	Plastic Bottle
5 Kg	Plastic Bottle
500 Gms	Plastic Bottle

Available Packages

500 Gms

1 Kg

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com