



64912

Diethylamine extrapure AR, 99.5%

Part A

CAS : 109-89-7

Molecular Formula : C₄H₁₁N

Molecular Weight : 73.14

Specifications

Appearance (Clarity)	Clear
Appearance (Form)	Liquid
Assay (GC)	min. 99.5%
Density (g/ml) @ 20°C	0.703-0.705
Refractive Index (20°C)	1.385-1.387
Boiling Range	55-57°C
Non Volatile Matter	max. 0.001%
Iron (Fe)	max. 0.00005%
Heavy Metals (Pb)	max. 0.0001%
Copper (Cu)	max. 0.00001%
Water (KF)	max. 0.1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Restricted for export from India
HSN Code	
2500 ml	29211200 (GST 18%)
500 ml	29211200 (GST 18%)
Type of Packing	
1000 ml	Glass Bottle
250 ml	Glass Bottle
2500 ml	Glass Bottle
500 ml	Glass Bottle

Available Packages

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com