



65163

Dextran Phenyl 40 (30-50)

Part B

Other Information

Application

(Phenyl dextrans may be used in the preparation of hydrophobic gels and also hydrophobic coatings.)

General Information

Storage

0 to 4°C (Freeze)

HSN Code

1 Gms

35051090 (GST 18%)

5 Gms

35051090 (GST 18%)

Available Packages

1 Gms

5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com