



65186

Albumin Bovine (pH 7) fraction V Reagent Grade (Bovine Serum Albumin, BSA), 99%

Part B

CAS : 9048-46-8

Specifications

Appearance (Colour)	Pale yellow to greenish yellow
Appearance (Form)	Crystalline powder
Solubility (Colour) 5% aq. solution	Pale yellow to greenish yellow
Solubility (Turbidity) 5% aq. solution	Clear
Assay (On dry basis)	min. 99%
pH (1% aq. solution)	6.5-7.5
Absorbance of 1% aq. Solution	
@280nm	6.34 - 6.60
Nitrogen (N)	min. 15.4%
Water (KF)	max. 5%
Fatty Acid	max. 0.05%
Protease	Not detected

Other Information

Application

Suitable for Binding & Transport Studies and all applications requiring low-fatty acid impurities. Also suitable for ELISA applications, Immunocytochemistry, Serology, as a Hapten Carrier and in Protein control systems.

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	38229090 (GST 12%)
100 Gms	38229090 (GST 12%)
25 Gms	38229090 (GST 12%)
Type of Packing	
10 Gms	Plastic Bottle
100 Gms	Plastic Bottle

Available Packages

10 Gms

25 Gms

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
