



65954

o-Tolidine extrapure AR, ExiPlus, Multi-Compendial, 98%

Part B

CAS : 119-93-7

Molecular Formula : C₁₄H₁₆N₂

Molecular Weight : 212.30

Specifications

ExiPlus	meets compendial specs of USP, BP, Ph.Eur
Appearance (Colour)	White to pale grey
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 2% solution in 5N HCl	Clear
Assay	min 98%
Melting Point	129 - 131°C
Loss on drying	max 0.2%
Sulphated Ash	max 0.1%

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	UN No.:3077 , IMCO Class No.:9 , Packing Group:III
HSN Code	
100 Gms	29215990 (GST 18%)
25 Gms	29215990 (GST 18%)
500 Gms	29215990 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

Available Packages

25 Gms

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com