



66120

Dimethyl Sulphoxide (DMSO) Dried, 99.5%, water 0.005%

Part A

CAS : 67-68-5

Molecular Formula : C₂H₆SO

Molecular Weight : 78.13

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 10
Assay (GC)	min. 99.5%
Density (g/ml) @ 20°C	1.099 - 1.101
Refractive Index (20°C)	1.478 - 1.479
Boiling Range	189 - 192°C
Freezing Point	17 - 18°C
Non Volatile Matter	max. 0.002%
Acidity (CH ₃ COOH)	max. 0.001%
Substance Discoloured by H ₂ SO ₄	Passes test
Water (KF)	max. 0.005%

Other Information

Notes

Packed with Dr. Dried Seal which are designed with corrosion free Crown Cap & with specially developed coloured septa. This enable repetitive needle puncturing easy and securing the consistency of the product from outside atmosphere.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
500 ml	29309099 (GST 18%)

Available Packages

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
