



66289

Aluminium Oxide Activated (Acidic) (Alumina)

Part A

CAS : 1344-28-1

Molecular Formula : Al2O3

Molecular Weight : 101.96

Specifications		
Appearance (Colour)	Almost white	
Appearance (Form)	Powder	
pH (10% aq. suspension)	~4.5	
Activity Grade (Brockmann)	1	
Sulphate (SO4)	max. 0.01%	
Iron (Fe)	max. 0.01%	
Heavy Metals (Pb)	max. 0.01%	
Adsorption Capacity	min. 1.0 mg o-Nitroaniline/gm	
Particle size (60 -325 Mesh BSS)	min. 90%	
Water Soluble Matter	max. 0.2%	

General Information		
Storage	25 to 40°C (Room Temperature)	
Shelf Life	60 Months	
IMDG Identification	Not Regulated for Transport (Non-Haz)	
HSN Code		
25 Kg	28182010 (GST 18%)	
5 Kg	28182010 (GST 18%)	
500 Gms	28182010 (GST 18%)	
Type of Packing		
5 Kg	Plastic Jar	
500 Gms	Plastic Bottle	

Available Packages	
500 Gms	
5 Kg	
25 Kg	

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com