



66400

Paraffin Wax w/ Ceresin Blocks, 58-60°C

Part A

CAS : 8002-74-2

Specifications

Appearance (Colour)	White to colourless mass, molten substance does not show fluorescence in daylight.
Appearance (Form)	Block
Odour	Odourless
Solubility (Turbidity) 10% solution in warm Xylene	Almost clear
Solubility (Turbidity) 10% solution in warm Chloroform	Almost clear
Acidity or alkalinity	Passes Test
Solidification point	58° - 60°C
Sulphated Ash	max 0.05%
Suitability for Histology	Passes Test

General Information

Storage	25 to 40°C (Room Temperature)
HSN Code	
500 Gms	27122090 (GST 18%)
1 Kg	27122090 (GST 18%)
2 Kg	27122090 (GST 18%)

Available Packages

500 Gms
1 Kg
2 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com