



66525

## 1,6-Diphenyl-1,3,5-Hexatriene (DPH) pure, 98%

Part B

CAS : 1720-32-7

Molecular Formula : C<sub>18</sub>H<sub>16</sub>

Molecular Weight : 232.32

### Specifications

Appearance (Colour)	Yellow
Appearance (Form)	Crystalline compound
Assay (GC)	min. 98.0%
Melting Point	199 - 203°C

### Other Information

Application

Used as an alternative to Nile red for fluorescent detection of lipid droplets in studies of membrane interiors and fluidity.

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
HSN Code	
1 Gms	29029090 (GST 18%)
5 Gms	29029090 (GST 18%)
Type of Packing	
5 Gms	Glass Bottle

### Available Packages

1 Gms

5 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com