



66532

Wang Resin, 1% DVB, 0.4-2mmol/g, 100-200 mesh

Part B

CAS : 201058-08-4

Other Information

Application

Suitable for peptide synthesis.
Commonly used linker in solid-phase peptide and small molecule synthesis.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	39039090 (GST 18%)
25 Gms	39039090 (GST 18%)
5 Gms	39039090 (GST 18%)

Available Packages

5 Gms
25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com