



66691

## Potassium Tellurite Anhydrous extrapure, 98%

Part A

CAS : 7790-58-1

Molecular Formula : K<sub>2</sub>TeO<sub>3</sub>

Molecular Weight : 253.80

### Specifications

Appearance (Colour)	White
Appearance (Form)	Powder
Solubility (Turbidity) 1% aq. solution	Clear
Solubility (Colour) 1% aq. solution	Colourless
Assay	min. 98%
Loss on drying	max. 1%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3284 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
10 Gms	28429010 (GST 18%)
100 Gms	28429010 (GST 18%)
25 Gms	28429010 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

### Available Packages

10 Gms  
25 Gms  
100 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com