



66886

## DEPC-Treated Water

Part E

CAS : 7732-18-5

### Specifications

Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Appearance (Clarity)	Clear
Absence of nucleases	Incubation of nucleic acids in DEPC treated water should show no degradation in comparison to the control as visualized after agarose gel electrophoresis.

### Other Information

Description	Water is deionized, diethylpyrocarbonate (DEPC) treated and 0.22 µm membrane-filtered. It is recommended to be used in any application where RNA is involved
-------------	--

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	38229090 (GST 12%)
500 ml	38229090 (GST 12%)

### Available Packages

100 ml  
500 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com