



67077

Ammonium Dichromate extrapure AR, ACS, 99.5%

Part A

CAS : 7789-09-5

Molecular Formula : $(\text{NH}_4)_2\text{Cr}_2\text{O}_7$

Molecular Weight : 252.06

Specifications

Appearance (Colour)	Orange - red
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Assay (T)	min. 99.5%
Insoluble matter	max. 0.003%
Loss on drying	max. 3%
Chloride (Cl)	max. 0.001%
Sulphate (SO ₄)	max. 0.002%
Iron (Fe)	max. 0.002%
Heavy Metals (Pb)	max. 0.005%
Calcium (Ca)	max. 0.002%
Sodium (Na)	max. 0.005%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1439 , IMCO Class No.:5.1 , Packing Group:II
HSN Code	
5 Kg	28415090 (GST 18%)
500 Gms	28415090 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

Available Packages

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com