



67209

Polysorbate 80 (Tween 80) for cell culture, Endotoxin (BET) 0.05EU/mg

Part B1

CAS : 9005-65-6

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Pale yellow to yellow to amber yellow
Appearance (Form)	viscous liquid
Density (g/ml) @ 20°C	1.077 - 1.081 g/ml
Refractive Index (20°C)	1.471 - 1.475
Acid value	max. 2.0
Hydroxyl value	65 - 80
Saponification value	45 - 55
Water (KF)	max. 3.0%
Suitability for Cell Culture	Passes test
Bacterial Endotoxin Test (BET) (EU/mg = IU/mg)	0.05 EU/mg

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
500 Gms	34024900 (GST 18%)

Available Packages

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com