



67258

## Polyguluronic Acid extrapure, 80%

Part B

CAS : 36562-70-6

Molecular Formula : C<sub>6</sub>H<sub>10</sub>O<sub>7</sub>

Molecular Weight : 194.14

### Specifications

Appearance (Colour)	White to Yellow to Light brown
Appearance (Form)	Powder
Assay	min. 80.0%
Water (KF)	max. 10.0%

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
HSN Code	
1 Gms	29389090 (GST 18%)
5 Gms	29389090 (GST 18%)

### Available Packages

1 Gms  
5 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com