



67299

## Zeatin mixed Isomers, ~95%

Part B

CAS : 13114-27-7

### Specifications

|                     |                      |
|---------------------|----------------------|
| Appearance (Colour) | White to off white   |
| Appearance (Form)   | Crystalline powder   |
| Assay               | ~95% (mixed isomers) |
| Melting Point       | 208 - 210°C          |

### General Information

|                     |                                       |
|---------------------|---------------------------------------|
| Storage             | 2 to 8°C (Refrigerate)                |
| Shelf Life          | 60 Months                             |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code            |                                       |
| 100 Mg              | 98020000 (GST 18%)                    |
| 25 Mg               | 98020000 (GST 18%)                    |
| Type of Packing     |                                       |
| 100 Mg              | Glass Bottle                          |

### Available Packages

25 Mg

---

100 Mg

---

250 Mg

---

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com