



67490

Tricyclohexylphosphine (Cy3P) extrapure, 94%

Part B

CAS : 2622-14-2

Molecular Formula : (C₆H₁₁)₃P

Molecular Weight : 280.43

Specifications

Appearance (Colour)	White
Appearance (Form)	Compound
Assay	min. 97%
Melting Point	81 - 83°C

Other Information

Application

Generally used in the following chemical reactions namely Heck Reaction, Hiyama, Sonogashira Coupling, Suzuki-Miyaura Coupling.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29314990 (GST 18%)
25 Gms	29314990 (GST 18%)
5 Gms	29314990 (GST 18%)

Available Packages

5 Gms
25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com