



67819

Ticarcillin Supplement (1mg/vl)

Part D

Specifications

Appearance (Colour)	White, fine powder
Solubility	Soluble in distilled water
Cultural Response	Inoculate and incubate at 25-30°C for 24-48 hours (Aseptically, add the rehydrated contents of 1 vial of Ticarcillin Supplement (67819) as per direction)
Organism	Yersinia enterocolitica ATCC 27729
Inoculum (cfu) 10-100	Growth : Good
Organism	Escherichia coli ATCC 25922
Inoculum	Growth : Inhibited
Organism	Staphylococcus aureus ATCC 25923
Inoculum	Growth : Inhibited

Other Information

Applications

Used for enrichment of Yersinia enterocolitica.

Composition

Ticarcillin	1.00mg/vial
-------------	-------------

Directions

1. Aseptically add 2 ml of sterile distilled water to 1 vial.
2. Mix well and add aseptically to sterile cooled media at 45-50°C and dispense as desired.
3. Add 2 ml (1 vial) of supplement to 988 ml of ITC Broth Base containing 10 ml of potassium chlorate supplement.

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 vl	38229090 (GST 12%)

Available Packages

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com