



67902

Tetraethylammonium Bromide (TEAB) pure, 98%

Part A

CAS : 71-91-0

Molecular Formula : C₈H₂₀NBr

Molecular Weight : 210.16

Specifications

Appearance (Colour)	White to pale yellow
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless to pale yellow
Assay	min. 98%
pH (10% aq. solution)	4.0 - 7.0

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kg	29239000 (GST 18%)
100 Gms	29239000 (GST 18%)
250 Gms	29239000 (GST 18%)
500 Gms	29239000 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

Available Packages

100 Gms

250 Gms

500 Gms

1 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com