



67911

Strontium Carbonate extrapure, 98%

Part A

CAS : 1633-05-2

Molecular Formula : SrCO₃

Molecular Weight : 147.63

Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Amorphous powder
Assay (NaOH titration)	min. 98%
Chloride (Cl)	max. 0.002%
Sulphate (SO ₄)	max. 0.1%
Iron (Fe)	max. 0.002%
Heavy Metals (Pb)	max. 0.002%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
250 Gms	28369200 (GST 18%)
500 Gms	28369200 (GST 18%)
Type of Packing	
250 Gms	Plastic Bottle
500 Gms	Plastic Bottle

Available Packages

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com