



67912

Cupric Acetate Monohydrate extrapure AR, 99%

Part A

CAS : 6046-93-1

Molecular Formula : (CH₃COO)₂Cu.H₂O

Molecular Weight : 199.65

Specifications

Appearance (Colour)	Dark green
Appearance (Form)	Crystalline compound
Solubility (Turbidity) 5% solution in hot 5N acetic acid	Clear
Solubility (Colour) 5% solution in hot 5N acetic acid	Blue
Assay	min. 99%
pH (5% aq. solution)	5.0 - 5.5
Chloride (Cl)	max. 0.001%
Sulphate (SO ₄)	max. 0.005%
Iron (Fe)	max. 0.002%
Heavy Metals (Pb)	max. 0.005%
Potassium (K)	max. 0.01%
Sodium (Na)	max. 0.01%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3077 , IMCO Class No.:9 , Packing Group:III
HSN Code	
250 Gms	29152990 (GST 18%)
Type of Packing	
250 Gms	Plastic Bottle

Available Packages

250 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com