



67959

Isooctane Fuel Grade for HPLC & UV Spectroscopy, 99.9%

Part A

CAS : 540-84-1

Molecular Weight : 816.64

Molecular Formula : C₃₅H₃₀Cl₄FeP₂Pd

General Information

Storage	25 to 40°C (Room Temperature)
HSN Code	
500 ml	29011000 (GST 18%)
2500 ml	29011000 (GST 18%)
25 Ltr	29011000 (GST 18%)

Available Packages

500 ml
2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com