



68030

## Lead (II) Nitrate ACS, 99%

Part A

CAS : 10099-74-8

Molecular Formula : Pb(NO<sub>3</sub>)<sub>2</sub>

Molecular Weight : 331.21

### Specifications

Appearance (Colour)	White
Appearance (Form)	Solid
Assay (EDTA Titration)	min. 99.0%
Insoluble matter	max. 0.005%
Chloride (Cl)	max. 0.001%
Calcium (Ca)	max. 0.005%
Copper (Cu)	max. 0.002%
Iron (Fe)	max. 0.001%
Potassium (K)	max. 0.005%
Sodium (Na)	max. 0.02%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1469 , IMCO Class No.:5.1(6.1) , Packing Group:II
HSN Code	
5 Kg	28342990 (GST 18%)
500 Gms	28342990 (GST 18%)

### Available Packages

500 Gms

5 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com