



68069

## SWCNT Type 10 Carbon Nanotubes Single-Walled in aq. media extrapure

Part C

CAS : 308068-56-6

### General Information

Storage	25 to 40°C (Room Temperature)
HSN Code	
20 ml	28030090 (GST 18%)

### Available Packages

20 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com