



68313

Phenol Red ACS, ExiPlus, Multi-Compendial

Part A

CAS : 143-74-8

Molecular Formula : C₁₉H₁₄O₅S

Molecular Weight : 354.38

Specifications

ExiPlus	meets compendial specs of BP, Ph.Eur
Appearance (Colour)	Dark red to Dark brown
Appearance (Form)	Powder
Loss on drying	max. 1%
Absorbance (A) of 1% solution (pH 8.0) in a 1cm cell v/s buffer @560nm	1100 - 1200
Clarity of solution	Passes test
Visual transition interval	Passes test

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	32049000 (GST 18%)
25 Gms	32049000 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

Available Packages

25 Gms

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com