



68434

## Mercuric Nitrate Monohydrate extrapure, 98%

Part A

CAS : 7783-34-8

Molecular Formula :  $\text{Hg}(\text{NO}_3)_2 \cdot \text{H}_2\text{O}$

Molecular Weight : 342.62

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% solution in 5N HNO <sub>3</sub>	Clear
Solubility (Colour) 10% solution in 5N HNO <sub>3</sub>	Colourless
Assay (T)	min. 98%
Chloride (Cl)	max. 0.005%
Sulphate (SO <sub>4</sub> )	max. 0.01%
Iron (Fe)	max. 0.002%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1625 , IMCO Class No.:6.1 , Packing Group:II
HSN Code	
100 Gms	28521000 (GST 18%)
25 Gms	28521000 (GST 18%)
250 Gms	28521000 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle
250 Gms	Glass Bottle

### Available Packages

25 Gms  
100 Gms  
250 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com