



68451

Deoxyribonucleic Acid Sodium Salt (DNA Sodium Salt Highly Polymerized) ex. Calf Thymus extrapure

Part B

CAS : 73049-39-5

Specifications

Appearance (Colour)	White fibrous
Appearance (Form)	Substance
Nitrogen (N)	~13.0 %
Phosphorus content	~8.0%
Absorbance (A) of 1% solution (pH 7.0) in a 1 cm cell v/s buffer @260nm	130 - 180
Electrophoresis test	Passes test

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Mg	29349990 (GST 18%)
25 Mg	29349990 (GST 18%)
250 Mg	29349990 (GST 18%)

Available Packages

25 Mg

100 Mg

250 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com