



68465

MWCNT Type 3 Carbon Nanotubes Multi Walled, Length 10-30µm, OD 10-20nm, 95%

Part C

CAS : 308068-56-6

Specifications

Appearance (Colour)	Black
Appearance (Form)	Powder
Assay	min. 95%
OD	10 - 20 nm
Length	10 - 30 µm

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1325 , IMCO Class No.:4.1 , Packing Group:III
HSN Code	
1 Gms	28030090 (GST 18%)
100 Gms	28030090 (GST 18%)
25 Gms	28030090 (GST 18%)
5 Gms	28030090 (GST 18%)
Type of Packing	
1 Gms	Plastic Bottle
100 Gms	Plastic Bottle
25 Gms	Plastic Bottle
5 Gms	Plastic Bottle

Available Packages

1 Gms
5 Gms
25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com