



68733

4-Fluorophenylboronic Acid extrapure, 98%

Part B

CAS : 1765-93-1

Molecular Formula : C₆H₆BF₂O₂

Molecular Weight : 139.92

Specifications

Appearance (Colour)	Off-white to pale yellow
Appearance (Form)	Crystalline powder
Assay (NT)	min. 98.0%
Melting range	260 - 265°C

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29420090 (GST 18%)
100 Gms	29420090 (GST 18%)
25 Gms	29420090 (GST 18%)
5 Gms	29420090 (GST 18%)

Available Packages

1 Gms
5 Gms
25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com