



68811

Diethylene Glycol (Digol) extrapure AR, 99%

Part A

CAS : 111-46-6

Molecular Formula : C₄H₁₀O₃

Molecular Weight : 106.12

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 10
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	1.114-1.116
Refractive Index (20°C)	1.446-1.447
Sulphated Ash	max. 0.01%
Iron (Fe)	max. 0.00005%
Heavy Metals (Pb)	max. 0.00005%
Ethylene Glycol	max. 0.75%
Water (KF)	max. 0.1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
2500 ml	29053990 (GST 18%)
500 ml	29053990 (GST 18%)
Type of Packing	
2500 ml	Glass Bottle
500 ml	Glass Bottle

Available Packages

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com