



68977

Rabbit anti-Bovine IgG (whole serum) (2mg/ml)

Part E

Specifications

Appearance (Form)	Liquid
Appearance (Colour)	Brown
Appearance (Clarity)	Clear
Activity	Tested in ELISA

Other Information

Description

Antiserum is produced in Rabbit using highly purified bovine IgG as the immunogen. Rabbit Anti-Bovine IgG is supplied as a liquid without any preservative. The antiserum is determined to be immunospecific for bovine IgG by immunoelectrophoresis(IEP) versus normal bovine serum and purified bovine IgG.

Application

ImmunoFluorescence (IF)
ImmunoCytoChemistry (ICC)
ImmunoHistoChemistry (IHC)
ELISA
Western Blotting
In Situ Hybridization

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 ml	38229090 (GST 12%)

Available Packages

1 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com