



68989

Trimethyl Orthoformate (TMOF) pure, 98%

Part B

CAS : 149-73-5

Molecular Formula : C₄H₁₀O₃

Molecular Weight : 106.12

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 98%
Density (g/ml) @ 20°C	0.968-0.969
Refractive Index (20°C)	1.378-1.379
Boiling Range	101-102°C

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3272 , IMCO Class No.:3 , Packing Group:II
HSN Code	
25 Ltr	29151300 (GST 18%)
2500 ml	29151300 (GST 18%)
500 ml	29151300 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

Available Packages

500 ml

2500 ml

25 Ltr

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com