



69031

## Phenol:Chloroform:Isoamyl Alcohol (25:24:1) pH 8.0 for molecular biology

Part A

CAS : 136112-00-0

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Pale yellow
Appearance (Form)	Liquid
pH	~ 8.0
DNase, RNase, Protease	Not detected
Stabilizer (8-Hydroxyquinoline)	Stabilized with 0.1% 8-Hydroxyquinoline

### Other Information

Notes	Suitable for the extraction of proteins from crude nucleic acid preparations
-------	--

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	UN No.:2810 , IMCO Class No.:6.1 , Packing Group:II
HSN Code	
100 ml	38229090 (GST 12%)
250 ml	38229090 (GST 12%)
500 ml	38229090 (GST 12%)

### Available Packages

100 ml  
250 ml  
500 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com