



69400

**Phospholipase D (Lecithinase D) ex.
Streptomyces chromofuscus, 40U/mg**

Part B

CAS : 9001-87-0

Molecular Weight : 57,000.00

Specifications

Appearance (Colour)	Yellowish brown
Appearance (Form)	Powder
Activity	min. 40 units/mg solid powder
Unit Definition	One unit will liberate 1.0 micromole of choline from L-?-phosphatidylcholine (egg yolk) per hour at pH 5.6 at 30°C.
Enzyme solution for activity determination	Prepared in Tris HCl buffer, pH 5.6.

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
200 Units	35079099 (GST 18%)
400 Units	35079099 (GST 18%)

Available Packages

200 Units

400 Units

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com