



69553

Eugenol extrapure, 99%

Part B

CAS : 97-53-0

Molecular Formula : C₁₀H₁₂O₂

Molecular Weight : 164.20

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	1.064 - 1.070 g/ml
Refractive Index (20°C)	1.540 - 1.542

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	29095090 (GST 18%)
500 ml	29095090 (GST 18%)

Available Packages

100 ml

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com