



69981

Benzene Dried, 99.7%, water 0.005%

Part A

CAS : 71-43-2

Molecular Formula : C6H6

Molecular Weight : 78.11

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max.10
Assay (GC)	min. 99.7%
Density (g/ml) @ 20°C	0.878-0.879
Refractive Index (20°C)	1.499-1.501
Boiling Range	80-81°C
Freezing Point	5°C
Non Volatile Matter	max. 0.001%
Acidity (HCl)	max. 0.001%
Substance Discoloured by H2SO4	Passes
Sulphur Compound (S)	max. 0.005%
Thiophene (C4H4S)	max. 0.0001%
Water (KF)	max. 0.005%

Other Information

Notes

Packed with Dr. Dried Seal which are designed with corrosion free Crown Cap & with specially developed coloured septa. This enable repetitive needle puncturing easy and securing the consistency of the product from outside atmosphere.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1114 , IMCO Class No.:3 , Packing Group:II
HSN Code	
1000 ml	29022000 (GST 18%)

Available Packages

1000 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com